KBC 11/15 - LAR2015

A - Aiguille 18G x 15cm	E - Mandrin perforé Ø2.3mm x 18.9cm
B - Stylet (mandrin)	F - Canule dentée Ø2.35mm x 24cm
C - Broche guide Ø0.9mm x 30cm	G - Poignée (amovible)
D - Gaine externe 11G x 15cm	H - Pousse carotte Ø1.7mm x 27.4cm

INDICATIONS

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PRECAUTIONS

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KBC 11/15 - LAR2015

A - Anaesthesia needle 18G x 15cm	E - Coaxial needle Ø2.3mm x 18.9cm
B - Stylet	F - Biopsy needle with sharp trephined tip Ø2.35mm x 24cm
C - Guiding stylet Ø0.9mm x 30cm	G - Removal plastic handle
D - Cannula 11G x 15cm	H - Ejector pin Ø1.7mm x 27.4cm

INDICATIONS FOR USE

Biopsy of bone lesions.

PRECAUTIONS

The device should be used only by or under supervision of a physician trained in biopsy or radiological interventions.

The device is supplied sterile. Check that the sterile package is intact. The device is single use and must not be cleaned and/or sterilized again. The incurred risk is a deterioration of the device and/or a contamination of the patient.

Carefully check the contents of the package prior to use and verify that all parts are present and undamaged.

Be careful not to bend the components of the device and use finger force only during the procedure. Never use a hammer.

When using other tools, first check the compatibility (length and gauge) of these tools with the device components.

After use, each component of the device may be a potential biohazard. Handle and dispose such devices in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

CONTRAINDICATIONS

The Coaxial Bone Biopsy Kit is contraindicated for use on patients who are receiving heavy anticoagulant therapy or who have a severe bleeding problem.

DIRECTIONS FOR USE

The following is a suggested method for penetration of the bone. The puncture towards a lesion should be accurately planned using diagnostic imaging.

After disinfection, apply local anesthetic using the anesthesia needle.

- Fig. 1** Perform anaesthesia with the anesthesia needle (A) all the way down to the periosteum or the lesion itself.
- Fig. 2** Insert the guiding stylet (C) in the anaesthesia needle (A) until bone is reached.
- Fig. 3** Hold firmly the guiding stylet (C) against the bone and remove the anaesthesia needle (A).
- Fig. 4** Place the Cannula and the Coaxial needle both assembled (D+E) over the Guiding stylet and insert it until bone is reached.
- Fig. 5** Hold firmly the Cannula (D) against the bone and remove the Guiding stylet (C) and the Coaxial needle (E).
- Fig. 6** Hold the Biopsy needle (F) and adjust the depth gauge (centimetric breakers) on the indicator wanted. The selected indicator will represent the protrusion of the Biopsy needle (F) beyond the tip of the Cannula (D).
- Fig. 7** Insert the Biopsy needle (F) together with its handle (G) in the Cannula (D). Carefully rotate the Biopsy needle clockwise until the depth gauge of the Biopsy needle reaches the Cannula Luer end. Then rotate counterclockwise and remove the Biopsy needle (F-E) from the Cannula. Remove the handle (G) from the Biopsy needle; insert the ejector pin (H) through the Biopsy needle (from the Luer end) and eject the sample.

Note 1: A syringe can be connected to the Luer end of the Biopsy needle and light aspiration may be applied during removal of the Biopsy needle.

Note 2: The Kensington Guide Wire Coaxial Bone Biopsy Kit is not intended for passage through thick cortical or hard bone. In such case use imperatively the Hard Bone Introducer (KTP11/15 - LAR3015).

STORAGE

Store in a dry, cool, dark place.
Do not use after expiration date.

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LAURANE Medical LLC