

# Kensington™



## GUIDE WIRE BONE BIOPSY KIT

This minimally invasive coaxial kit features:

- ▶ Seldinger placement of introducer via 20G guide wire
- ▶ Fewer scans and minimal trauma to soft tissue
- ▶ Ideal for deep bone or narrow access procedures
- ▶ Proven high diagnostic and sampling success rates<sup>1</sup>
- ▶ Superior lytic and sclerotic sampling
- ▶ Compatible ingress for supplemental ablation
- ▶ Safe access for disc procedures with 11G needle

### Contents:

- anesthesia needle
- guide wire
- hollow stylet/tapered-tip cannula introducer
- trephine-tip biopsy needle
- ejector pin



Kensington reference codes:

Provided sterile, one unit per box

Code	Reference	Gauge	Length	Core
KBC 11/10	Kensington 11/10	11 gauge	10 cm	1.8 mm
KBC 11/15	Kensington 11/15	11 gauge	15 cm	1.8 mm



[www.lauranemedical.com](http://www.lauranemedical.com)

<sup>1</sup>Laredo J-D, Hamze B, Jeribi R. Percutaneous biopsy of osteoid osteomas prior to percutaneous treatment using two different biopsy needles. Cardiovasc Intervent Radiol. 2009 Sep;32(5):998-1003. doi: 10.1007/s00270-009-9635-2