

# Madison™

## COMPREHENSIVE BONE BIOPSY KIT


This all-inclusive coaxial kit features:

- ▶ Unparalleled manual hard bone access
- ▶ Proven high diagnostic and sampling success rates<sup>1</sup>
- ▶ Laurane's unique perforating introducer cannula
- ▶ Direct access and easy handling
- ▶ High-grade steel for outstanding durability
- ▶ Compatible ingress for supplemental ablation

Contents (in 2 trays):

- trocar-tip stylet/perforating-tip cannula introducer
- exchangeable drill insert
- trephine-tip biopsy needle
- ejector pin

An animated video is available on our website: [www.lauranemedical.com](http://www.lauranemedical.com)



Madison reference codes:

Provided sterile, one unit per box

Code	Reference	Gauge	Length	Core
KDPD 11/6.5	Madison 11/6.5	11 gauge	6.5 cm	1.8 mm
KDPD 11/8.5	Madison 11/8.5	11 gauge	8.5 cm	1.8 mm
KDPD 11/10	Madison 11/10	11 gauge	10 cm	1.8 mm
KDPD 11/15	Madison 11/15	11 gauge	15 cm	1.8 mm
KDPD 13/9*	Madison 13/9	13 gauge	9 cm	1.4 mm
KDP 13/6*†	Madison Mini 13/6	13 gauge	6 cm	1.4 mm

\*Madison 13 gauge kits feature a removable, screw-lock handle

†Madison Mini is a single tray and does not include a drill



[www.lauranemedical.com](http://www.lauranemedical.com)