

# Westbrook™

## BASIC BONE BIOPSY KIT

This cost-effective coaxial kit features:

- ▶ Direct access and comfortable handling
- ▶ Fast, efficient penetration of normal density bone
- ▶ Proven high diagnostic and sampling success rates<sup>1</sup>
- ▶ Superior sclerotic and lytic sampling
- ▶ Compatible ingress for ablation techniques
- ▶ High-grade steel for outstanding durability
- ▶ Trephine-tip biopsy needle minimising crush artefacts

### Contents:

- trocar-tip stylet/tapered-tip cannula introducer
- trephine-tip biopsy needle
- ejector pin



Westbrook reference codes:

Provided sterile, one unit per box

Code	Reference	Gauge	Length	Core
KBD 11/10	Westbrook 11/10	11 gauge	10 cm	1.8 mm
KBD 11/15	Westbrook 11/15	11 gauge	15 cm	1.8 mm



[www.lauranemedical.com](http://www.lauranemedical.com)

<sup>1</sup>Laredo J-D, Hamze B, Jeribi R. Percutaneous biopsy of osteoid osteomas prior to percutaneous treatment using two different biopsy needles. Cardiovasc Intervent Radiol. 2009 Sep;32(5):998-1003. doi: 10.1007/s00270-009-9635-2